



painaway
Home Care Laser

Support Positive
Recovery for
Workplace and
Auto Accident
Injuries

The only FDA cleared super pulsed laser therapy for home use.

The cordless, portable, effective pain relief device that is research proven, safe, and easy to use. The PainAway Home Care Laser provides temporary pain relief and enhances post-operative care for injured workers' comp or personal injury patients, all in the comfort of their homes.

“ As a doctor, I’ve been treating patients for 13 years. In our offices, we’ve been using laser technology to treat acute & chronic pain. I highly recommend the laser because it will help your pain go away. ”

Nosson S. Goldfarb, MD | Integrative Medicine



Easy to Use

Simple interface and treatment modes



Opioid-Free Pain Relief

Effective for acute and chronic pain



FDA OTC Cleared

For patient use at home



Convenient

Portable / cordless, 8 hour rechargeable battery



Safe to Use

No side effects and quick 5-minute treatment times



Proof is in the Research

A collection of summaries and scientific papers on safety and effectiveness of Multi Radiance Medical's Super Pulsed Laser Therapy in different fields of application.



FDA-cleared alternative to pain medications

- ✓ Low level laser therapy has been found to offer superior pain-relieving effects compared to other electrotherapeutic modalities in both chronic injuries and the early stages of acute injuries.
- ✓ Multi Radiance Medical technology creates a high photon density, strongly reducing pain and improving microcirculation.
- ✓ Portable device treats muscle, tendon, ligament, connective tissue, bone, nerve, and dermal tissues in a non-invasive, drug-free modality.
- ✓ The super pulsed laser produces high powered light in billionth-of-a-second pulses that drives photons deep into the target tissue.

PainAway Home Care Laser Program

With our full-service program we offer 3rd-party billing services for payer for pre-authorization. Once the payer pre-authorizes the workmans' comp or personal injury referral, we will contact the physician's clinic/practice to inform them of the approval and ship the PainAway Home Care device to the clinic/practice to dispense to the patient, or we can send the device directly to the patient's home. If the referral is denied due to missing information or other reasons, we will contact the clinic/practice to obtain the information for re-submission.



ELBOW PAIN

NECK PAIN

ANKLE PAIN



SHOULDER PAIN

WRIST PAIN

LEG PAIN



Accelerate Rehabilitation

Treatment protocols will have your patients returning to their activities of daily living faster and in less pain, reducing the need for prescription pain medications.

Easy setup, easy operation.

The PainAway Home Care Laser is safe, effective, and easy to use. It is cordless with a rechargeable battery for treatments at home or on the go.

Mode	Treatment
1	For pain immediately following an acute injury used up to 48 hours post-injury.
2	For pain exceeding 48 hours in duration. Ideal for muscles and tendons.
3	Provides maximum pain relief for most injuries.
All treatments last for 5 minutes	

The PainAway Home Care device has three (3) treatment modes to choose from, providing 15W of Peak Super Pulsed Power with pulsing Red and Infrared Radiances. Patients can treat one or two times per day.



Want to improve outcomes, satisfaction, and your bottom line? Please contact your rep or visit multiradiance.com.



multiradiance.com
(877) 542.0761



Info@VastlandSupply.com
877-VASTLAND / 877-827-8526