



The only FDA-cleared super pulsed laser therapy for home use.

The cordless, portable, and effective pain relief device that is research proven, safe, and easy to use. The PainAway Laser 2.0 can provide temporary pain relief and enhance post-operative care for patients, all in the comfort of their homes.

“ As a doctor, I’ve been treating patients for 13 years. In our offices, we’ve been using laser technology to treat acute & chronic pain. I highly recommend the laser because it will help your pain go away. ”

Nosson S. Goldfarb, MD | Integrative Medicine



Easy to Use

Simple interface and 5 treatment modes



Opioid-Free Pain Relief

Effective for acute and chronic pain with no side effects



FDA OTC Cleared

For patient use at home



Convenient

Portable / cordless with rechargeable battery



Fast Results

Quick treatments in 5 minutes or less



TECHNOLOGY THAT HAS BEEN VALIDATED IN VITRO, IN VIVO, IN CONTROLLED CLINICAL TRIALS AND IN CLINICAL PRACTICE.

Positive and prompt recoveries are an important step towards getting back to the things that matter. Multi Radiance laser technology has been used safely and effectively used among a wide range of patient populations for years. Physical therapists, orthopedics, outpatient clinics, and chiropractic offices, as well as professional and collegiate teams worldwide use Multi Radiance laser therapy every day.

FDA-cleared alternative to pain medications

- ✓ Low level laser therapy has been found to offer superior pain-relieving effects compared to other electrotherapeutic modalities in both chronic injuries and the early stages of acute injuries.
- ✓ Multi Radiance Medical technology creates a high photon density, strongly reducing pain and improving microcirculation.
- ✓ Non-invasive, drug-free modality that treats muscles, tendons, ligaments, connective tissues, bones, nerves, and dermal tissues.
- ✓ The super pulsed laser produces high powered light in billionth-of-a-second pulses that drives photons deep into the target tissue.

Patients benefit most when they take on a significant role in optimal management of their own pain.



Using the Device

The PainAway device is simple to use, has five (5) treatment modes to choose from, and provides 25W of Peak Super Pulsed Power along with pulsing Red and Infrared radiances. Patients can treat one or two times per day per physician recommendation, with treatment sessions lasting 1-5 minutes.





Accelerate Rehabilitation

Treatment protocols will have your patients returning to their activities of daily living faster and in less pain, reducing the need for prescription pain medications.



LUMBAR PAIN



KNEE PAIN



LEG PAIN



NECK PAIN

Proof is in the Research

A collection of summaries and scientific papers on safety and effectiveness of Multi Radiance Medical's Super Pulsed Laser Therapy in different fields of application.



SHOULDER PAIN



WRIST PAIN



ANKLE PAIN



ELBOW PAIN

painaway 2.0 LASER 2.0


Support your patient's injury recovery, rehabilitation, and pain management, even when they're not in your office.



Easy setup, easy operation.

The PainAway Laser 2.0 is safe, effective, and easy to use. It is cordless with a rechargeable battery for treatments at home or on the go.

Mode	Treatment	Time
1	Increasing Circulation, Inflammation	1 min
2	Joint Pain, Arthritis	1 min
3	Pain and Relaxing Muscles	3 mins
4	Acute Injury (Under 48 hrs only)	2 mins
5	General and Chronic Pain	5 mins



The PainAway Laser 2.0 device has five (5) treatment modes to choose from, providing 25W of Peak Super Pulsed Power with pulsing Laser, Red Light and Infrared Radiances. Patients can treat one or two times per day.

Want to improve outcomes, satisfaction, and your bottom line? Please contact your rep or visit multiradiance.com.